

Model	SinoAngio-I200
Syringe	Single syringe
Size of syringe	150 ML sterile syringe
Injection Volume	0.1ML - Volume of syringe
Injection Delay	0-600 s
Injection Phases	1-8 Phases
Injection Protocol	Upto 120 protocols
Tilt Lockout	Rotate head to arm-down, minimize risk of air embolish
Control Console	12.1 inch real color touch LCD for both local and remote (Remote is optional)
ISI	Synchronization of injection delay and scanning dealy
Pressure limit	100 - 1200 Psi
Flow rate	0.1mL/s`50.0mL/s step: 0.1mLs/ for 0.1mL/ssXs10.0mL/s 1mL/s for 10mL/s X<50.0mL/s